


Commitment to Standard

When you need a sheath that follows accepted clinical practice for vessel puncture or incision and guidewire insertion, you need the performance of **BIOTEQ[®]** sheath introducer set.

- Meet, even exceed ISO 11070 standard “Sterile single use intravascular catheter introducer” in all aspects.

Easily Accessible

Smooth transition
between dilator
and sheath



Color coded sheath
hub and dilator,
with rotatable
suture eye



Easy-to-lock
dilator and
sheath hub




Hub with durable
hemostasis valve
inside



Sideport with
stopcock for infusion
and contrast
injection



Wide variety of sets and sheaths



DESCRIPTION AND FEATURES

- Wide variety of sets and sheaths in adult and pediatric sizes
- Smooth transition between sheath and dilator
- Maximum flexibility without kinking or compression
- Taper tip design to reduce potential tissue trauma
- Lubricious sheath and dilator
- Radiopaque sheath with hemostasis valve and sideport with attached 3-way stopcock

ORDERING INFORMATION:

Kit Type	Model No.	Inner Diameter	Sheath Length	Compatible Guidewire	Guidewire Type	Puncture Needle	Syringe	Scalpel	Color Code
A Kit in pouch without guidewire	BT-SI-1-040110-03150-00	4 Fr	11 cm	0.035"	-	-	-	-	Red
	BT-SI-1-050110-03150-00	5 Fr	11 cm	0.038"	-	-	-	-	Gray
	BT-SI-1-060110-03150-00	6 Fr	11 cm	0.038"	-	-	-	-	Green
	BT-SI-1-070110-03150-00	7 Fr	11 cm	0.038"	-	-	-	-	Orange
	BT-SI-1-040230-03150-00	4 Fr	23 cm	0.035"	-	-	-	-	Red
	BT-SI-1-050230-03150-00	5 Fr	23 cm	0.038"	-	-	-	-	Gray
	BT-SI-1-060230-03150-00	6 Fr	23 cm	0.038"	-	-	-	-	Green
B Kit in pouch with guidewire	BT-SI-1-040110-43150-00	4 Fr	11 cm	0.035"	45 cm J	-	-	-	Red
	BT-SI-1-050110-53150-00	5 Fr	11 cm	0.038"	45 cm J	-	-	-	Gray
	BT-SI-1-060110-53150-00	6 Fr	11 cm	0.038"	45 cm J	-	-	-	Green
	BT-SI-1-070110-53150-00	7 Fr	11 cm	0.038"	45 cm J	-	-	-	Orange
	BT-SI-1-040230-43150-00	4 Fr	23 cm	0.035"	80 cm J	-	-	-	Red
	BT-SI-1-050230-53150-00	5 Fr	23 cm	0.038"	80 cm J	-	-	-	Gray
	BT-SI-1-060230-53150-00	6 Fr	23 cm	0.038"	80 cm J	-	-	-	Green
C Kit in pouch with guidewire & needle	BT-SI-1-040110-43150-66	4 Fr	11 cm	0.035"	45 cm J	18 G X 3"	-	-	Red
	BT-SI-1-050110-53150-66	5 Fr	11 cm	0.038"	45 cm J	18 G X 3"	-	-	Gray
	BT-SI-1-060110-53150-66	6 Fr	11 cm	0.038"	45 cm J	18 G X 3"	-	-	Green
	BT-SI-1-070110-53150-66	7 Fr	11 cm	0.038"	45 cm J	18 G X 3"	-	-	Orange
	BT-SI-1-040230-43150-66	4 Fr	23 cm	0.035"	80 cm J	18 G X 3"	-	-	Red
	BT-SI-1-050230-53150-66	5 Fr	23 cm	0.038"	80 cm J	18 G X 3"	-	-	Gray
	BT-SI-1-060230-53150-66	6 Fr	23 cm	0.038"	80 cm J	18 G X 3"	-	-	Green
D Kit in tray with guidewire, needle & syringe	BT-SI-1-040110-43150-33	4 Fr	11 cm	0.035"	45 cm J	18 G X 3"	3 c.c.	-	Red
	BT-SI-1-050110-53150-33	5 Fr	11 cm	0.038"	45 cm J	18 G X 3"	3 c.c.	-	Gray
	BT-SI-1-060110-53150-33	6 Fr	11 cm	0.038"	45 cm J	18 G X 3"	3 c.c.	-	Green
	BT-SI-1-070110-53150-33	7 Fr	11 cm	0.038"	45 cm J	18 G X 3"	3 c.c.	-	Orange
	BT-SI-1-040230-43150-33	4 Fr	23 cm	0.035"	80 cm J	18 G X 3"	3 c.c.	-	Red
	BT-SI-1-050230-53150-33	5 Fr	23 cm	0.038"	80 cm J	18 G X 3"	3 c.c.	-	Gray
	BT-SI-1-060230-53150-33	6 Fr	23 cm	0.038"	80 cm J	18 G X 3"	3 c.c.	-	Green
E Kit in tray with guidewire, needle, syringe & scalpel	BT-SI-1-040110-43150-44	4 Fr	11 cm	0.035"	45 cm J	18 G X 3"	3 c.c.	Yes	Red
	BT-SI-1-050110-53150-44	5 Fr	11 cm	0.038"	45 cm J	18 G X 3"	3 c.c.	Yes	Gray
	BT-SI-1-060110-53150-44	6 Fr	11 cm	0.038"	45 cm J	18 G X 3"	3 c.c.	Yes	Green
	BT-SI-1-070110-53150-44	7 Fr	11 cm	0.038"	45 cm J	18 G X 3"	3 c.c.	Yes	Orange
	BT-SI-1-040230-43150-44	4 Fr	23 cm	0.035"	80 cm J	18 G X 3"	3 c.c.	Yes	Red
	BT-SI-1-050230-53150-44	5 Fr	23 cm	0.038"	80 cm J	18 G X 3"	3 c.c.	Yes	Gray
	BT-SI-1-060230-53150-44	6 Fr	23 cm	0.038"	80 cm J	18 G X 3"	3 c.c.	Yes	Green
BT-SI-1-070230-53150-44	7 Fr	23 cm	0.038"	80 cm J	18 G X 3"	3 c.c.	Yes	Orange	



For more information, please contact:
sales@bioteq.com.tw

BIOTEQ CORPORATION
 5F-6, No.23, Sec.1, Chang-An East Rd.,
 Taipei, 104, Taiwan
 TEL:886-2-25710269
 FAX:886-2-25361967

ISO 13485, GMP, CE Certified
 ©Bioteq Corp. July, 2012
 Brochure Number:BT- SI-01